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

	Page
Erection	9
1.1. Conditions to be imposed on the Place of Erection	9
1.2. Lighting	9
1.3. Adjustment of the Horizontal Position	9
1.4. Loosening the Transport Safety Device	9
1.5. Special Features	9
2. Adjusting the Optical Indication	10
2.1. Adjusting the Optical Indication	10
2.2. Zero Adjustment	10
2.3. Sensitivity checking and correcting	10
Operator Attention	11
3.1. Weighing without Tare Compensation	11
3.2. Weighing with the Aid of a Tare Compensation	11
3.3. Weighing with an Optical Micrometer	11
3.4. Weighing below the Balance	12
Service and Maintenance	12

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## 1. Erection

The balance is dispatched in a shock- and impact-resistant package. The balance and its accessories have to be taken out cautiously. Have regard to the marking "Top" !!!

### 1.1. Conditions to be imposed on the Place of Erection

- Vibrationless erection
- No intensive draught
- Maximum temperature differential within some hours  $\leq 5^\circ\text{C}$
- Prevention of any direct thermal influence

### 1.2. Lighting

- Take off the lid 4 (fig. 2).
- Set the voltage selector 2 (fig. 2) in accordance with the existing line voltage (set to 220 V by the manufacturer).
- Switch in the slide switch 3 (fig. 2) checking the functioning of the incandescent lamp.
- If the lamp does not shine, it is necessary to check the line voltage, the incandescent lamp, and the safety fuse in the voltage selector (safety fuse to be exchanged only at the "0"-position of the selector).

### 1.3. Adjustment of the Horizontal Position

- Place the balance at its definite location.
- Adjust the horizontal position by means of the front adjusting legs 7 (fig. 1) according to the air levels 8 (fig. 2).
- Lock the adjusting legs by means of the lock nuts 6 (fig. 1).


### 1.4. Loosening the Transport Safety Device

- Put the pan 2 (fig. 1) on. (In case of type 34.005 and 34.011 the upper part of the protective ring 5 (fig. 2) has to be slipped over).
- Unlock the balance by turning the hand wheel (below the inscription "Transport safety device") as far as the stop.
- Turn the operating knob 9 (fig. 1) until the zero will appear in the field 4 (fig. 1).

### 1.5. Special Features

In case of type 34.008 there have now to be taken into consideration the following particularities:

After loosening the transport safety device the enclosed pointer 3 (fig. 4) is screwed on to the end of the lever, the balance pan being depressed. The complete top structure 2 (fig. 4) is put on vertically from above and fastened by means of the clamping screws 4 (fig. 4).

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The scale image is looked at in the reading mirror (1) (fig. 4) which is adjustable in height.

In finding out the weight value one has to take into consideration the added load indicated in field 4 (fig. 1).  
 On the right-hand border of the reading field there is provided a coarse survey of the inclination range (0...100 g). Readjust pointer 3 (fig. 4) only when it is absolutely necessary.  
 The reading mirror shall be cleaned only by means of a soft-brush or cloth.

**2. Adjustment**

**2.1. Adjusting the Optical Indication**

- After switching on the incandescent lamp by means of the slide switch 3 (fig. 2), effect the critical focusing of the scale image by turning the knob 7 (fig. 2).
  - Only if required, adjust the maximum brightness of the scale image in the reading field 5 (fig. 1) by adjusting the lamp holder 1 (fig. 2).
- If the balance does not swing through, turn the operating knob 9 (fig. 1) through.

**2.2. Zero Adjustment**

Before adjusting or checking the zero point there are always necessary the following manipulations:

- Bring the indexing weights to the zero position by means of operating knob 9 (fig. 1).
- Set the tare compensating device 11 (fig. 1) to "0" as far as the perceptible stop according to the survey indication 3 (fig. 1).
- Set the micrometer to "0" by means of operating knob 6 (fig. 2).

Operate the adjusting wheel 3 (fig. 3) or 5 (fig. 4) until the graduation line of the scale marked by "0" will coincide with the reading mark.

The front adjusting legs can also be employed for the zero point regulation. In doing that, pay attention to the air level (the bubble must remain within the marking).

If the zero point cannot be longer adjusted in that way, then a correction on the lever will be necessary:

- Remove the balance pan 2 (fig. 1) and the pan support 1 (fig. 1) (lift the pan support off, in case of type 34.010 screw it off).
- Loosen the screws 10 (fig. 1).
- Take off the casing 8 (fig. 1). In case of type 34.008 it is necessary in the meantime to remove the top structure 2 (fig. 4).
- Put on again the pan support and the pan.

**2.3. Sensitivity checking and correcting**

After every zero point correction by means of the taring nut, the sensitivity has to be checked (see the following survey):

- Load the balance with the weight X.
- The graduation line of the scale marked by Y must coincide with the reading mark.
- Permissible deviation: Z

Type	X	Y	Z
34.001 and 34.002	100 g	100	± 0,1 g
34.003 and 34.004	10 g	10	± 0,01 g
34.008	1000 g	1000	± 1,0 g
34.005 and 34.011	1000 g	1000	± 1,0 g
34.010	1000 g	1000	± 1,0 g


In case of an insufficient deflection set the adjusting nut 2 (fig. 3) higher, and in case of an extensive deflection set it lower.

**3. Operator Attention**

- Before the weighing operation check the zero point (zero-point checking)
- Place the material to be weighed cautiously in order to spare the knife-edges and knife-edge supports.

**3.1. Weighing without Tare Compensation**

The weight value of a material to be weighed which is within the inclination range can be immediately read off on the ground-glass plate 5 (fig. 1) or on the reading mirror 1 (fig. 4). In case of greater weight values the operating knob 9 (fig. 1) has to be turned on until the scale will balance out.

Type 34.001

**3.2. Weighing with the Aid of a Tare Compensation**

In case of balances with built-in tare compensating device masses can be tared out steplessly. The survey indication 3 (fig. 1) will permit a coarse survey.


By taking the pan 2 (fig. 1) off, a certain mass can be additionally tared without operating the indexing mechanism.

In the following survey there are represented the tare compensating possibilities of the various balance types.

Type	Stepless tare compensation	Weight of the pan	Max. tare compensation
34.002	0...100 g	200 g	300 g
34.004	0...100 g	100 g	200 g
34.008	0...100 g	500 g	1000 g
34.011	0...500 g	500 g	1000 g
34.010	0...500 g	500 g	1000 g

Performing the weighing operation with tare compensation:

- Tare compensation of receptacles put on by means of operating knob 11 (fig. 1) until the zero point will be reached
- to be continued as under 3.1.


Type 34.002

**3.3. Weighing with an Optical Micrometer**


- As under 3.1. or 3.2. respectively.
- By turning the operating knob 6 (fig. 2) make the next graduation line of the scale coincide with the reading mark.

Unambiguous numerical display of the weighing result in the reading field.

Type 34.004

### 3.4. Weighing below the Balance

- After loosening the 2 screws swing the lid on the bottom side of the balance away. In case of a possible tilting of the balance, the transport safety device has to be set to "locked".
- Attach the material to be weighed to the hook provided for that purpose.
- Tare compensation by taking off the balance pan or by operating the tare compensating device.

### 4. Servicing and Maintenance

When transporting the balance, one must by all means take care that first of all the indexing weights are locked by means of the operating knob 9 (fig. 1), which is shown by the display of a red arrow in the viewing window 4 (fig. 1). After that the transport safety device has to be set to "locked". Transport safety device to be operated only in case of an important change of site, not between the successive weighing operations.

In order to increase the service life of the incandescent lamps, it is advisable to switch the lighting off during weighing intervals of some length.

Though, as a matter of principle, the balance does not require any maintenance owing to its technical functionality, we recommend to have an inspection or readjustment performed approximately every 2 years.

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